

Antibodies Product Datasheet

Product Name:	IFN- α 1 Antibody
Catalog Number:	C30079
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IFNA1/13. AA range:140-189
Concentration:	1 mg/ml
Gene Name:	IFNA1/IFNA13
Protein Name:	Interferon alpha-1/13
Other Name:	IFNA1; IFNA13; Interferon alpha-1/13; IFN-alpha-1/13; Interferon alpha-D; LeIF D
Human Swiss Prot No.:	P01562
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFN- α 1 Polyclonal Antibody detects endogenous levels of IFN- α 1 protein.
Storage and Stability:	-20°C/1 year