

Antibodies Product Datasheet

Product Name:	CGR19 Antibody
Catalog Number:	C30085
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human CGRRF1. AA range:283-332
Concentration:	1 mg/ml
Gene Name:	CGRRF1
Protein Name:	Cell growth regulator with RING finger domain protein 1
Other Name:	CGRRF1; CGR19; RNF197; Cell growth regulator with RING finger domain protein 1; Cell growth regulatory gene 19 protein; RING finger protein 197
Human Swiss Prot No.:	Q99675
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CGR19 Polyclonal Antibody detects endogenous levels of CGR19 protein.
Storage and Stability:	-20°C/1 year