

Antibodies Product Datasheet

Product Name:	INDOL1 Antibody
Catalog Number:	C30096
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IDO2. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	IDO2
Protein Name:	Indoleamine 2,3-dioxygenase 2
Other Name:	IDO2; INDOL1; Indoleamine 2,3-dioxygenase 2; IDO-2; Indoleamine 2,3-dioxygenase-like protein 1; Indoleamine-pyrrole 2,3-dioxygenase-like protein 1
Human Swiss Prot No.:	Q6ZQW0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	INDOL1 Polyclonal Antibody detects endogenous levels of INDOL1 protein.
Storage and Stability:	-20°C/1 year