

Antibodies Product Datasheet

Product Name:	MK Antibody
Catalog Number:	C30098
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human MDK. AA range:94-143
Concentration:	1 mg/ml
Gene Name:	MDK
Protein Name:	Midkine
Other Name:	MDK; MK1; NEGF2; Midkine; MK; Amphiregulin-associated protein; ARAP; Midgestation and kidney protein; Neurite outgrowth-promoting factor 2; Neurite outgrowth-promoting protein
Human Swiss Prot No.:	P21741
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MK Polyclonal Antibody detects endogenous levels of MK protein.
Storage and Stability:	-20°C/1 year