

Antibodies Product Datasheet

Product Name:	Lamin B1 Antibody
Catalog Number:	C30101
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LMNB1. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	LMNB1
Protein Name:	Lamin-B1
Other Name:	LMNB1; LMN2; LMNB; Lamin-B1
Human Swiss Prot No.:	P20700
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Lamin B1 Polyclonal Antibody detects endogenous levels of Lamin B1 protein.
Storage and Stability:	-20°C/1 year