

Antibodies Product Datasheet

Product Name:	IL-12A p35 Antibody
Catalog Number:	C30103
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL12A. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	IL12A
Protein Name:	Interleukin-12 subunit alpha
Other Name:	IL12A; NKSF1; Interleukin-12 subunit alpha; IL-12A; Cytotoxic lymphocyte maturation factor 35 kDa subunit; CLMF p35; IL-12 subunit p35; NK cell stimulatory factor chain 1; NKSF1
Human Swiss Prot No.:	P29459
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-12A p35 Polyclonal Antibody detects endogenous levels of IL-12A p35 protein.
Storage and Stability:	-20°C/1 year