

Antibodies Product Datasheet

Product Name:	Angptl1/2 Antibody
Catalog Number:	C30112
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ANGPTL1/2. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	ANGPTL1/ANGPTL2
Protein Name:	Angiopoietin-related protein 1/Angiopoietin-related protein 2
Other Name:	ANGPTL1; ANG3; ANGPT3; ARP1; PSEC0154; Angiopoietin-related protein 1; Angiopoietin-3; ANG-3; Angiopoietin-like protein 1; ANGPTL2; ARP2; Angiopoietin-related protein 2; Angiopoietin-like protein 2
Human Swiss Prot No.:	O95841
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Angptl1/2 Polyclonal Antibody detects endogenous levels of Angptl1/2 protein.
Storage and Stability:	-20°C/1 year