

Antibodies Product Datasheet

Product Name:	CUL-4B Antibody
Catalog Number:	C30113
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CUL4B. AA range:711-760
Concentration:	1 mg/ml
Gene Name:	CUL4B
Protein Name:	Cullin-4B
Other Name:	CUL4B; KIAA0695; Cullin-4B; CUL-4B
Human Swiss Prot No.:	Q13620
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CUL-4B Polyclonal Antibody detects endogenous levels of CUL-4B protein.
Storage and Stability:	-20°C/1 year