

Antibodies Product Datasheet

Product Name:	Cdc34 Antibody
Catalog Number:	C30115
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CDC34. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	CDC34
Protein Name:	Ubiquitin-conjugating enzyme E2 R1
Other Name:	CDC34; UBCH3; UBE2R1; Ubiquitin-conjugating enzyme E2 R1; Ubiquitin-conjugating enzyme E2-32 kDa complementing; Ubiquitin-conjugating enzyme E2-CDC34; Ubiquitin-protein ligase R1
Human Swiss Prot No.:	P49427
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc34 Polyclonal Antibody detects endogenous levels of Cdc34 protein.
Storage and Stability:	-20°C/1 year