

Antibodies Product Datasheet

Product Name:	FGF-5 Antibody
Catalog Number:	C30117
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FGF5. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	FGF5
Protein Name:	Fibroblast growth factor 5
Other Name:	FGF5; Fibroblast growth factor 5; FGF-5; Heparin-binding growth factor 5; HBGF-5; Smag-82
Human Swiss Prot No.:	P12034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-5 Polyclonal Antibody detects endogenous levels of FGF-5 protein.
Storage and Stability:	-20°C/1 year