

Antibodies Product Datasheet

Product Name:	SCYL1BP1 Antibody
Catalog Number:	C30128
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human GORAB. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GORAB
Protein Name:	RAB6-interacting golgin
Other Name:	GORAB; NTKLBP1; SCYL1BP1; RAB6-interacting golgin; N-terminal kinase-like-binding protein 1; NTKL-BP1; NTKL-binding protein 1; hNTKL-BP1; SCY1-like 1-binding protein 1; SCYL1-BP1; SCYL1-binding protein 1
Human Swiss Prot No.:	Q5T7V8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SCYL1BP1 Polyclonal Antibody detects endogenous levels of SCYL1BP1 protein.
Storage and Stability:	-20°C/1 year