

Antibodies Product Datasheet

Product Name:	Stra8 Antibody
Catalog Number:	C30132
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human STRA8. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	STRA8
Protein Name:	Stimulated by retinoic acid gene 8 protein homolog
Other Name:	STRA8; Stimulated by retinoic acid gene 8 protein homolog
Human Swiss Prot No.:	Q7Z7C7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Stra8 Polyclonal Antibody detects endogenous levels of Stra8 protein.
Storage and Stability:	-20°C/1 year