

Antibodies Product Datasheet

Product Name:	ASIC3 Antibody
Catalog Number:	C30145
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ASIC3. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	ASIC3
Protein Name:	Acid-sensing ion channel 3
Other Name:	ASIC3; ACCN3; SLNAC1; TNAC1; Acid-sensing ion channel 3; ASIC3; hASIC3; Amiloride-sensitive cation channel 3; Neuronal amiloride-sensitive cation channel 3; Testis sodium channel 1; hTNaC1
Human Swiss Prot No.:	Q9UHC3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ASIC3 Polyclonal Antibody detects endogenous levels of ASIC3 protein.
Storage and Stability:	-20°C/1 year