

Antibodies Product Datasheet

Product Name:	Annexin VII Antibody
Catalog Number:	C30151
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ANXA7. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	ANXA7
Protein Name:	Annexin A7
Other Name:	ANXA7; ANX7; SNX; OK/SW-cl.95; Annexin A7; Annexin VII; Annexin-7; Synexin
Human Swiss Prot No.:	P20073
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Annexin VII Polyclonal Antibody detects endogenous levels of Annexin VII protein.
Storage and Stability:	-20°C/1 year