

Antibodies Product Datasheet

Product Name:	CEACAM1/5 Antibody
Catalog Number:	C30155
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CEACAM1/CEACAM5. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CEACAM1/CEACAM5
Protein Name:	Carcinoembryonic antigen-related cell adhesion molecule 1/Carcinoembryonic antigen-related cell adhesion molecule 5
Other Name:	CEACAM1; BGP; BGP1; Carcinoembryonic antigen-related cell adhesion molecule 1; Biliary glycoprotein 1; BGP-1; CD66a; CEACAM5; CEA; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; CEA; Meconium antigen 100; CD66e
Human Swiss Prot No.:	P13688
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CEACAM1/5 Polyclonal Antibody detects endogenous levels of CEACAM1/5 protein.
Storage and Stability:	-20°C/1 year