

Antibodies Product Datasheet

Product Name:	ECA39 Antibody
Catalog Number:	C30158
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BCAT1. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	BCAT1
Protein Name:	Branched-chain-amino-acid aminotransferase, cytosolic
Other Name:	BCAT1; BCT1; ECA39; Branched-chain-amino-acid aminotransferase, cytosolic; BCAT(c); Protein ECA39
Human Swiss Prot No.:	P54687
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ECA39 Polyclonal Antibody detects endogenous levels of ECA39 protein.
Storage and Stability:	-20°C/1 year