

Antibodies Product Datasheet

Product Name:	E-Selectin Antibody
Catalog Number:	C30167
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300.WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human SELE. AA range:100-150
Concentration:	1 mg/ml
Gene Name:	SELE ELAM1
Protein Name:	E-Selectin
Other Name:	selectin E
Human Swiss Prot No.:	P16581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	E-Selectin Polyclonal Antibody detects endogenous levels of E-Selectin
Storage and Stability:	-20°C/1 year