

Antibodies Product Datasheet

Product Name:	MIP-1 β Antibody
Catalog Number:	C30171
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CCL4. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	CCL4
Protein Name:	C-C motif chemokine 4
Other Name:	CCL4; LAG1; MIP1B; SCYA4; C-C motif chemokine 4; G-26 T-lymphocyte-secreted protein; HC21; Lymphocyte activation gene 1 protein; LAG-1; MIP-1-beta(1-69);Macrophage inflammatory protein 1-beta; MIP-1-beta; PAT 744; Protein H400; SIS-gamma; Small-inducible P13236
Human Swiss Prot No.:	P13236
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MIP-1 β Polyclonal Antibody detects endogenous levels of MIP-1 β protein.
Storage and Stability:	-20°C/1 year