

Antibodies Product Datasheet

Product Name:	BLC Antibody
Catalog Number:	C30177
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CXCL13. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	CXCL13
Protein Name:	C-X-C motif chemokine 13
Other Name:	CXCL13; BCA1; BLC; SCYB13; C-X-C motif chemokine 13; Angie; B cell-attracting chemokine 1; BCA-1; B lymphocyte chemoattractant; CXC chemokine BLC; Small-inducible cytokine B13
Human Swiss Prot No.:	O43927
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BLC Polyclonal Antibody detects endogenous levels of BLC protein.
Storage and Stability:	-20°C/1 year