

Antibodies Product Datasheet

Product Name:	Dkk-1 Antibody
Catalog Number:	C30180
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human DKK1. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	DKK1
Protein Name:	Dickkopf-related protein 1
Other Name:	DKK1; Dickkopf-related protein 1; Dickkopf-1; Dkk-1; hDkk-1; SK
Human Swiss Prot No.:	O94907
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dkk-1 Polyclonal Antibody detects endogenous levels of Dkk-1 protein.
Storage and Stability:	-20°C/1 year