

Antibodies Product Datasheet

Product Name:	LIF Antibody
Catalog Number:	C30217
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human LIF. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	LIF
Protein Name:	Leukemia inhibitory factor
Other Name:	LIF; HILDA; Leukemia inhibitory factor; LIF; Differentiation-stimulating factor; D factor; Melanoma-derived LPL inhibitor; MLPLI; Emfilermin
Human Swiss Prot No.:	P15018
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LIF Polyclonal Antibody detects endogenous levels of LIF protein.
Storage and Stability:	-20°C/1 year