

Antibodies Product Datasheet

Product Name:	D3DR Antibody
Catalog Number:	C30226
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from D3DR . at AA range: 181-230
Concentration:	1 mg/ml
Gene Name:	DRD3
Protein Name:	D(3) dopamine receptor
Other Name:	DRD3; D(3) dopamine receptor; Dopamine D3 receptor
Human Swiss Prot No.:	P35462
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	D3DR Polyclonal Antibody detects endogenous levels of D3DR protein.
Storage and Stability:	-20°C/1 year