

Antibodies Product Datasheet

Product Name:	IL-5R α Antibody
Catalog Number:	C30234
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human IL5RA. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	IL5RA
Protein Name:	Interleukin-5 receptor subunit alpha
Other Name:	IL5RA; IL5R; Interleukin-5 receptor subunit alpha; IL-5 receptor subunit alpha; IL-5R subunit alpha; IL-5R-alpha; IL-5RA; CDw125; CD125
Human Swiss Prot No.:	Q01344
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-5R α Polyclonal Antibody detects endogenous levels of IL-5R α protein.
Storage and Stability:	-20°C/1 year