

Antibodies Product Datasheet

Product Name:	KLK1 Antibody
Catalog Number:	C30238
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human KLK1. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	KLK1
Protein Name:	Kallikrein-1
Other Name:	KLK1; Kallikrein-1; Kidney/pancreas/salivary gland kallikrein; Tissue kallikrein
Human Swiss Prot No.:	P06870
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KLK1 Polyclonal Antibody detects endogenous levels of KLK1 protein.
Storage and Stability:	-20°C/1 year