

Antibodies Product Datasheet

Product Name:	MCM3 Antibody
Catalog Number:	C30244
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human MCM3. AA range:641-690
Concentration:	1 mg/ml
Gene Name:	MCM3
Protein Name:	DNA replication licensing factor MCM3
Other Name:	MCM3; DNA replication licensing factor MCM3; DNA polymerase alpha holoenzyme-associated protein P1; P1-MCM3; RLF subunit beta; p102
Human Swiss Prot No.:	P25205
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MCM3 Polyclonal Antibody detects endogenous levels of MCM3 protein.
Storage and Stability:	-20°C/1 year