

Antibodies Product Datasheet

Product Name:	Lipocalin-1 Antibody
Catalog Number:	C30247
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LCN1. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	LCN1
Protein Name:	Lipocalin-1
Other Name:	LCN1; VEGP; Lipocalin-1; Tear lipocalin; Tlc; Tear prealbumin; TP; Von Ebner gland protein; VEG protein
Human Swiss Prot No.:	P31025
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Lipocalin-1 Polyclonal Antibody detects endogenous levels of Lipocalin-1 protein.
Storage and Stability:	-20°C/1 year