

Antibodies Product Datasheet

Product Name:	NFATc1 Antibody
Catalog Number:	C30249
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human NFATC1. AA range:881-930
Concentration:	1 mg/ml
Gene Name:	NFATC1
Protein Name:	Nuclear factor of activated T-cells cytoplasmic 1
Other Name:	NFATC1; NFAT2; NFATC; Nuclear factor of activated T-cells, cytoplasmic 1; NF-ATc1; NFATc1; NFAT transcription complex cytosolic component; NF-ATc; NFATc
Human Swiss Prot No.:	O95644
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NFATc1 Polyclonal Antibody detects endogenous levels of NFATc1 protein.
Storage and Stability:	-20°C/1 year