

Antibodies Product Datasheet

Product Name:	GABA T-1 Antibody
Catalog Number:	C30272
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000, IHC-p 1:500-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SLC6A1. AA range:170-220
Concentration:	1 mg/ml
Gene Name:	SLC6A1 GABATR GABT1 GAT1
Protein Name:	Sodium- and chloride-dependent GABA transporter 1 (GAT-1) (Solute carrier family 6 member 1)
Other Name:	Sodium- and chloride-dependent GABA transporter 1 (GAT-1) (Solute carrier family 6 member 1)
Human Swiss Prot No.:	P30531
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous GABA T-1
Storage and Stability:	-20°C/1 year