

Antibodies Product Datasheet

Product Name:	DCI Antibody
Catalog Number:	C30286
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DCI . at AA range: c-terminal
Concentration:	1 mg/ml
Gene Name:	ECI1
Protein Name:	Enoyl-CoA delta isomerase 1 mitochondrial
Other Name:	ECI1; DCI; Enoyl-CoA delta isomerase 1, mitochondrial; 3,2-trans-enoil-CoA isomerase; Delta(3),Delta(2)-enoil-CoA isomerase; D3,D2-enoil-CoA isomerase; Dodecenoyl-CoA isomerase
Human Swiss Prot No.:	P42126
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DCI Polyclonal Antibody detects endogenous levels of DCI protein.
Storage and Stability:	-20°C/1 year