

Antibodies Product Datasheet

Product Name:	FMO3 Antibody
Catalog Number:	C30287
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FMO3. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	FMO3
Protein Name:	Dimethylaniline monooxygenase [N-oxide-forming] 3
Other Name:	FMO3; Dimethylaniline monooxygenase [N-oxide-forming] 3; Dimethylaniline oxidase 3; FMO II; FMO form 2; Hepatic flavin-containing monooxygenase 3; FMO 3; Trimethylamine monooxygenase
Human Swiss Prot No.:	P31513
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FMO3 Polyclonal Antibody detects endogenous levels of FMO3 protein.
Storage and Stability:	-20°C/1 year