

Antibodies Product Datasheet

Product Name:	HEXA Antibody
Catalog Number:	C30300
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from HEXA . at AA range: 121-170
Concentration:	1 mg/ml
Gene Name:	HEXA
Protein Name:	Beta-hexosaminidase subunit alpha
Other Name:	HEXA; Beta-hexosaminidase subunit alpha; Beta-N-acetylhexosaminidase subunit alpha; Hexosaminidase subunit A; N-acetyl-beta-glucosaminidase subunit alpha
Human Swiss Prot No.:	P06865
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HEXA Polyclonal Antibody detects endogenous levels of HEXA protein.
Storage and Stability:	-20°C/1 year