

Antibodies Product Datasheet

Product Name:	HLA-DM β Antibody
Catalog Number:	C30301
Reactivity:	Human
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 40-100
Concentration:	1 mg/ml
Gene Name:	HLA-DMB DMB RING7
Protein Name:	HLA class II histocompatibility antigen, DM beta chain (MHC class II antigen DMB) (Really interesting new gene 7 protein)
Other Name:	HLA class II histocompatibility antigen, DM beta chain (MHC class II antigen DMB) (Really interesting new gene 7 protein)
Human Swiss Prot No.:	P28068
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HLA-DM β
Storage and Stability:	-20°C/1 year