

Antibodies Product Datasheet

Product Name:	LASP-1 Antibody
Catalog Number:	C30302
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LASP1. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	LASP1
Protein Name:	LIM and SH3 domain protein 1
Other Name:	LASP1; MLN50; LIM and SH3 domain protein 1; LASP-1; Metastatic lymph node gene 50 protein; MLN 50
Human Swiss Prot No.:	Q14847
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LASP-1 Polyclonal Antibody detects endogenous levels of LASP-1 protein.
Storage and Stability:	-20°C/1 year