

Antibodies Product Datasheet

Product Name:	AACT Antibody
Catalog Number:	C30313
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SERPINA3. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	SERPINA3
Protein Name:	Alpha-1-antichymotrypsin
Other Name:	SERPINA3; AACT; GIG24; GIG25; Alpha-1-antichymotrypsin; ACT; Cell growth-inhibiting gene 24/25 protein; Serpin A3
Human Swiss Prot No.:	P01011
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AACT Polyclonal Antibody detects endogenous levels of AACT protein.
Storage and Stability:	-20°C/1 year