

Antibodies Product Datasheet

Product Name:	CD1A/B Antibody
Catalog Number:	C30314
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 220-270
Concentration:	1 mg/ml
Gene Name:	CD1A CD1B
Protein Name:	T-cell surface glycoprotein CD1A/B (CD antigen CD1A/B)
Other Name:	
Human Swiss Prot No.:	P06126/P29016
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous CD1A/B
Storage and Stability:	-20°C/1 year