

Antibodies Product Datasheet

Product Name:	CD1D Antibody
Catalog Number:	C30317
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Antigen-presenting glycoprotein CD1d at AA range: 161-210
Concentration:	1 mg/ml
Gene Name:	CD1D
Protein Name:	Antigen-presenting glycoprotein CD1d
Other Name:	CD1D; Antigen-presenting glycoprotein CD1d; R3G1; CD1d
Human Swiss Prot No.:	P15813
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD1D Polyclonal Antibody detects endogenous levels of CD1D protein.
Storage and Stability:	-20°C/1 year