

Antibodies Product Datasheet

Product Name:	CD2 Antibody
Catalog Number:	C30319
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CD2. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CD2
Protein Name:	T-cell surface antigen CD2
Other Name:	CD2; SRBC; T-cell surface antigen CD2; Erythrocyte receptor; LFA-2; LFA-3 receptor; Rosette receptor; T-cell surface antigen T11/Leu-5; CD2
Human Swiss Prot No.:	P06729
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD2 Polyclonal Antibody detects endogenous levels of CD2 protein.
Storage and Stability:	-20°C/1 year