

Antibodies Product Datasheet

Product Name:	CD3- δ Antibody
Catalog Number:	C30320
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD3D. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	CD3D
Protein Name:	T-cell surface glycoprotein CD3 delta chain
Other Name:	CD3D; T3D; T-cell surface glycoprotein CD3 delta chain; T-cell receptor T3 delta chain; CD3d
Human Swiss Prot No.:	P04234
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD3- δ Polyclonal Antibody detects endogenous levels of CD3- δ protein.
Storage and Stability:	-20°C/1 year