

Antibodies Product Datasheet

Product Name:	Integrin α L Antibody
Catalog Number:	C30325
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ITGAL. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	ITGAL
Protein Name:	Integrin alpha-L
Other Name:	ITGAL; CD11A; Integrin alpha-L; CD11 antigen-like family member A; Leukocyte adhesion glycoprotein LFA-1 alpha chain; LFA-1A; Leukocyte function-associated molecule 1 alpha chain; CD11a
Human Swiss Prot No.:	P20701
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin α L Polyclonal Antibody detects endogenous levels of Integrin α L protein.
Storage and Stability:	-20°C/1 year