

Antibodies Product Datasheet

Product Name:	Integrin α D Antibody
Catalog Number:	C30328
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ITGAD. AA range:901-950
Concentration:	1 mg/ml
Gene Name:	ITGAD
Protein Name:	Integrin alpha-D
Other Name:	ITGAD; Integrin alpha-D; ADB2; CD11 antigen-like family member D; Leukointegrin alpha D; CD11d
Human Swiss Prot No.:	Q13349
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin α D Polyclonal Antibody detects endogenous levels of Integrin α D protein.
Storage and Stability:	-20°C/1 year