

Antibodies Product Datasheet

Product Name:	CD21 Antibody
Catalog Number:	C30335
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CR2. AA range:381-430
Concentration:	1 mg/ml
Gene Name:	CR2
Protein Name:	Complement receptor type 2
Other Name:	CR2; C3DR; Complement receptor type 2; Cr2; Complement C3d receptor; Epstein-Barr virus receptor; EBV receptor; CD21
Human Swiss Prot No.:	P20023
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD21 Polyclonal Antibody detects endogenous levels of CD21 protein.
Storage and Stability:	-20°C/1 year