

Antibodies Product Datasheet

Product Name:	CD26 Antibody
Catalog Number:	C30338
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human DPP4. AA range:350-400
Concentration:	1 mg/ml
Gene Name:	DPP4
Protein Name:	Dipeptidyl peptidase 4
Other Name:	DPP4; ADCP2; CD26; Dipeptidyl peptidase 4; ADABP; Adenosine deaminase complexing protein 2; ADCP-2; Dipeptidyl peptidase IV; DPP IV; T-cell activation antigen CD26; TP103; CD26
Human Swiss Prot No.:	P27487
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD26 Polyclonal Antibody detects endogenous levels of CD26 protein.
Storage and Stability:	-20°C/1 year