

Antibodies Product Datasheet

Product Name:	CD32-A Antibody
Catalog Number:	C30342
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FCGR2A. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	FCGR2A
Protein Name:	Low affinity immunoglobulin gamma Fc region receptor II-a
Other Name:	FCGR2A; CD32; FCG2; FCGR2A1; IGFR2; Low affinity immunoglobulin gamma Fc region receptor II-a; IgG Fc receptor II-a; CDw32; Fc-gamma RII-a; Fc-gamma-RIIa; FcRII-a; CD32
Human Swiss Prot No.:	P12318
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD32-A Polyclonal Antibody detects endogenous levels of CD32-A protein.
Storage and Stability:	-20°C/1 year