

Antibodies Product Datasheet

Product Name:	CD42d Antibody
Catalog Number:	C30354
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human GP5. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	GP5
Protein Name:	Platelet glycoprotein V
Other Name:	GP5; Platelet glycoprotein V; GPV; Glycoprotein 5; CD42d
Human Swiss Prot No.:	P40197
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD42d Polyclonal Antibody detects endogenous levels of CD42d protein.
Storage and Stability:	-20°C/1 year