

Antibodies Product Datasheet

Product Name:	Integrin α 6 Antibody
Catalog Number:	C30361
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ITGA6. AA range:901-950
Concentration:	1 mg/ml
Gene Name:	ITGA6
Protein Name:	Integrin alpha-6
Other Name:	ITGA6; Integrin alpha-6; CD49 antigen-like family member F; VLA-6; CD49f
Human Swiss Prot No.:	P23229
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin α 6 Polyclonal Antibody detects endogenous levels of Integrin α 6 protein.
Storage and Stability:	-20°C/1 year