

Antibodies Product Datasheet

Product Name:	CD55 Antibody
Catalog Number:	C30366
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD55. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	CD55
Protein Name:	Complement decay-accelerating factor
Other Name:	CD55; CR; DAF; Complement decay-accelerating factor; CD55
Human Swiss Prot No.:	P08174
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD55 Polyclonal Antibody detects endogenous levels of CD55 protein.
Storage and Stability:	-20°C/1 year