

Antibodies Product Datasheet

Product Name:	CD58 Antibody
Catalog Number:	C30368
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD58. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CD58
Protein Name:	Lymphocyte function-associated antigen 3
Other Name:	CD58; LFA3; Lymphocyte function-associated antigen 3; Ag3; Surface glycoprotein LFA-3; CD58
Human Swiss Prot No.:	P19256
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD58 Polyclonal Antibody detects endogenous levels of CD58 protein.
Storage and Stability:	-20°C/1 year