

Antibodies Product Datasheet

Product Name:	CD66a/b/c Antibody
Catalog Number:	C30375
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 130-190
Concentration:	1 mg/ml
Gene Name:	CEACAM1/CEACAM6/CEACAM8
Protein Name:	Carcinoembryonic antigen-related cell adhesion molecule 1/6/8 (CD antigen CD66a/b/c)
Other Name:	
Human Swiss Prot No.:	P13688/P40199/P31997
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous CD66a/b/c
Storage and Stability:	-20°C/1 year