

Antibodies Product Datasheet

Product Name:	CD69 Antibody
Catalog Number:	C30379
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD69. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	CD69
Protein Name:	Early activation antigen CD69
Other Name:	CD69; CLEC2C; Early activation antigen CD69; Activation inducer molecule; AIM; BL-AC/P26; C-type lectin domain family 2 member C; EA1; Early T-cell activation antigen p60; GP32/28; Leukocyte surface antigen Leu-23; MLR-3; CD69
Human Swiss Prot No.:	Q07108
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD69 Polyclonal Antibody detects endogenous levels of CD69 protein.
Storage and Stability:	-20°C/1 year