

## Antibodies Product Datasheet

|                               |                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CD79b Antibody                                                                                                                                                          |
| <b>Catalog Number:</b>        | C30384                                                                                                                                                                  |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                       |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                                                  |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.                                                                    |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                        |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                       |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                   |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from the Internal region of human CD79B. AA range:61-110                                                 |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                 |
| <b>Gene Name:</b>             | CD79B                                                                                                                                                                   |
| <b>Protein Name:</b>          | B-cell antigen receptor complex-associated protein beta chain                                                                                                           |
| <b>Other Name:</b>            | CD79B; B29; IGB; B-cell antigen receptor complex-associated protein beta chain; B-cell-specific glycoprotein B29; Ig-beta; Immunoglobulin-associated B29 protein; CD79b |
| <b>Human Swiss Prot No.:</b>  | P40259                                                                                                                                                                  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                 |
| <b>Specificity:</b>           | CD79b Polyclonal Antibody detects endogenous levels of CD79b protein.                                                                                                   |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                            |